

PATIENT CARE SERIES – ICCU & ICU BEDS

ICU - MEM

- **ICCU BED ELECTRIC**
- 8 Function Bed
- CRCA Bed Top
- ABS Railings with inbuilt Membrane Switch.
- ABS Head Board with inbuilt Membrane Switch.
- Dual Wheel Castors
- Sturdy CRCA Construction.



ICU - EXS

- **ICCU BED ELECTRIC**
- 5 Function Bed
- CRCA Bed Top
- ABS Railings
- ABS Head Board with inbuilt Safety Ball Bumpers
- Dual Wheel Castors
- Sturdy CRCA Construction.



ICU - M

- **ICU Manual**
- 5 Function Bed
- CRCA Perforated Bed Top.
- ABS Railings (2Big + 2Small)
- ABS Head Board
- 125mm Wheel Castors
- Sturdy CRCA Construction.



Options / Add-ons

- **Wheels :** Double Wheel Castors or Central Locking Double Wheel Castors.
- **Bed Top :** 4 Sectional ABS Bed Top & X-Ray Cassette Holder in backrest.
 - **Nurse Control Panel.**
 - **Battery Backup.**
 - **Four Fold Mattress.**



Due to continuous development and innovations, we reserve the rights to change the design and look of the product without prior intimation.

PATIENT CARE BEDS SERIES – RECOVERY BEDS



RECOVERY-EX

- **RECOVERY BED ELECTRIC**
- 3 Function Bed
- Height/Back & Leg.
- ABS Bed Top
- ABS Railings (4Small)
- ABS Head Board
- Dual Wheel Castors
- Sturdy CRCA Construction.



RECOVERY-M6

- **RECOVERY BED MANUAL**
- 3 Function Bed
- Height/Back & Leg.
- CRCA Perforated Bed Top
- ABS Railings (4Small)
- ABS Head Board
- 125mm Wheel Castors
- Sturdy CRCA Construction.

Options / Add-ons

- **Wheels :** Double Wheel Castors or Central Locking Double Wheel Castors.
- **Bed Top :** 4 Sectional ABS Bed Top & X-Ray Cassette Holder in backrest.
 - Nurse Control Panel.
 - Battery Backup.
 - Four Fold Mattress.



Due to continuous development and innovations, we reserve the rights to change the design and look of the product without prior intimation.

PATIENT CARE BEDS SERIES – ROOM CARE BEDS

PATIENT CARE-02

- **FULL FOWLER BED ELECTRIC**
- 2 Function Bed
- Back & Leg.
- CRCA Bed Top
- Collapsible SS Side Railings / ABS Railings
- Wooden Head Board
- Sturdy CRCA Construction.
- Choice of Wheels



PATIENT CARE-11

- **FULL FOWLER BED MANUAL**
- 2 Function Bed
- Back & Leg.
- CRCA Bed Top
- Collapsible SS Side Railings.
- ABS Head Board
- Sturdy CRCA Construction.
- Choice of Wheels



COMFORT - 06

- **SEMI FOWLER BED MANUAL**
- Single Function Bed
- Back.
- CRCA Bed Top
- Collapsible SS Side Railings.
- ABS Head Board
- Sturdy CRCA Construction.
- Choice of Wheels



Options / Add-ons

- **Wheels** : Double Wheel Castors or Central Locking Double Wheel Castors.
- **ABS Railings or SS Collapsible Railings**
- **Mattress**



Due to continuous development and innovations, we reserve the rights to change the design and look of the product without prior intimation.