

MULTI PURPOSE SURGICAL TABLES

ULTIMA – ENDO-1

- **SURGICAL TABLE**
ENDOSCOPY-ELECTRIC
- Sturdy Stainless Steel Construction.
- Electrically operated Height, Tilting (Both Sides), Trend. (Both Sides) & Top Slide.
- Mattress Fixed Table Top.
- Provided with Standard Stainless Steel Accessories



ULTIMA – EXS

- **SURGICAL TABLE**
MULTISECTIONAL-ELECTRIC
- Sturdy Stainless Steel Construction.
- Electrically operated Height, Tilting (Both Sides), Trend. (Both Sides), Flex/Reflex & Top Slide.
- X-RAY Translucent Mattress Table Top.
- Divided Leg with Rotateable and droppable option.
- Manually Operated Kidney Bridge.
- Provided with Standard Stainless Steel Accessories



Due to continuous development and innovations, we reserve the rights to change the design and look of the product without prior intimation.