

MEDICINE CARTS/ TROLLIES

MEDCART CR4

- **MEDICINE CART**
- Sturdy CRCA Construction.
- Provided with 5 Drawers and 1 Open Shelf
- Stainless Steel Top Shelf.
- Dual Wheel Castors for extra stability.
- Iv Rod
- Corner Buffers



MEDCART

- **MEDICINE CART**
- Sturdy CRCA & Aluminium Alloy Frame Construction.
- Specially Designed Heavy Duty Moulded Cart Top with provision to hold IV Rod and Monitor Table
- Fully Covered with ABS Plastic Cladding with CRCA Double Support.
- 6 Medicine Drawers (3 Small and 3 Big size, to be run on extra smooth telescopic Channels.
- Central Locking system for Drawers
- Dual Wheel Castors for extra stability.
- Iv Rod
- Corner Buffers



CR-MED-05

- **CRASH CART TROLLEY**
- Sturdy Stainless Steel or CRCA Construction as desired.
- 2 Stainless Steel Shelves.
- 6 Plastic Medicine Drawers
- 6 Coloured Bins.
- CPR Table
- IV Rod
- Cylinder Cage



Due to continuous development and innovations, we reserve the rights to change the design and look of the product without prior intimation.