

DELIVERY BED SERIES



LDR EXT

- **LDR - DELIVERY BED ELECTRIC**
- 3 Function Delivery Bed
- Electric Height, Back and Trendeleburg.
- Fitted with X-Ray Translucent Bed Top.
- ABS Side Rails.
- Sliding Leg Section to convert Labour Table to Bed.
- Sturdy CRCA Construction.
- Provided with SS. Bowl Basin
- Fully Covered with FRP Moulded Covers. Supplied with SS Accessories Trolley and standard accessories



LDR -05

- **DELIVERY BED ELECTRIC**
- 3 Function Delivery Bed
- Electric Height, Back and Trendeleburg.
- Fitted with Fixed Mattress Bed Top
- SS Side Rails.
- Sturdy CRCA Construction.
- Provided with SS. Bowl Basin
- Base Covered with FRP Moulded Covers.



Due to continuous development and innovations, we reserve the rights to change the design and look of the product without prior intimation.

DELIVERY BED SERIES

GYNOX -11

- **DELIVERY BED**
- 3 Function Delivery Bed
- Electrically Operated Height, Back & Trendelenburg.
- CRCA Bed Top
- Collapsible ABS Side Railings.
- ABS Head Board
- Sturdy CRCA Construction.
- Provided with SS. Bowl Basin
- Optional : Beckalite or SS Bed Top.



GYNOX -02

- **DELIVERY BED**
- 3 Function Delivery Bed
- Electrically Operated Height, Backrest adjustable on Ratchet and Trendelenburg adjustable by Crank Handle
- CRCA Bed Top
- Collapsible SS Side Railings.
- ABS Head Board
- Sturdy CRCA Construction.
- Provided with SS. Bowl Basin
- Optional : Beckalite or SS Bed Top.
- Optional : ABS Railings



GYNOX – H6

- **DELIVERY BED**
- 3 Function Delivery Bed
- Foot Pump Operated Height, Backrest adjustable on Ratchet and Trendelenburg adjustable by Crank Handle
- CRCA Bed Top
- Collapsible SS Side Railings.
- ABS Head Board
- Sturdy CRCA Construction.
- Provided with SS. Bowl Basin
- Optional : Beckalite or SS Bed Top.
- Optional : ABS Railings



Due to continuous development and innovations, we reserve the rights to change the design and look of the product without prior intimation.