

# EXAMINATION TABLES AND COUCHES



## EXAMO-RGH

- **OPD EXAMINATION TABLE (Electric)**
- 2 Function Examination Table
- Electrically Operated Height Backrest adjustment.
- Optional : Gas Spring assisted Backrest also available instead of Electric.
- Fully Covered Mountings.
- 2 Sectional Fixed Mattress Top
- Sturdy CRCA Construction.
- Both side sliding Drawer for Instruments.



## EXAMO-PR

- **PREMIUM OPD EXAMINATION TABLE**
- Fixed Height OPD Examination Table.
- Option of Electrically Operated or Gas Spring assisted Backrest
- 2 Sectional Fixed Mattress Top
- Provided with Sliding Leg Section.
- Provided with Sliding Foot Stool.
- Provided with 5 Wooden Drawers.
- Sturdy CRCA Construction.



## EXAMO-OPD

- **OPD EXAMINATION TABLE**
- Fixed Height OPD Examination Table.
- Option of Electrically Operated or Gas Spring assisted Backrest.
- 2 Sectional Fixed Mattress Top
- Provided with 3 Wooden Drawers.
- Provided with Wooden Sliding Cabinet.
- Sturdy CRCA Construction.
- BP Apparatus Tray.



Due to continuous development and innovations, we reserve the rights to change the design and look of the product without prior intimation.

# EXAMINATION TABLES AND COUCHES



## EXAMO-GLT

- **ULTRASOUND EXAMINATION TABLE (Electric)**
- 2 Function Mobile Examination Table
- Electrically Operated Height & Backrest adjustment.
- Fixed Mattress Top
- Special 22" Min and 36" Max Height for Ultrasound comfort
- Base Back FRP Cover
- Sturdy CRCA Construction.

## EXAMO-UT

- **UROLOGY EXAMINATION TABLE (Electric)**
- 3 Function Mobile Examination Table
- Electrically Operated Height, Backrest & Trend. adjustment.
- Fixed Mattress Top
- Fitted with S.S. Lithotomy Crutches with PU Moulded Pads
- PU Arm rests.
- Provided with FRP Moulded Basin
- Base Back FRP Cover
- Sturdy CRCA Construction.

## EXAMO-GN2

- **GYNAE EXAMINATION TABLE Mobile (Electric)**
- 3 Function Mobile Examination Table
- Electrically Operated Height, Backrest & Legrest adjustment.
- Fixed Mattress Top
- Fitted with S.S. Lithotomy Crutches with PU Moulded Pads
- Provided with Stainless Steel Collapsible Railings.
- Sturdy CRCA Construction.
- **Optional** : Foot Switch for Electric Functions

