

PATIENT CARE BEDS SERIES – RECOVERY BEDS



RECOVERY-EX

- **RECOVERY BED ELECTRIC**
- 3 Function Bed
- Height/Back & Leg.
- ABS Bed Top
- ABS Railings (4Small)
- ABS Head Board
- Dual Wheel Castors
- Sturdy CRCA Construction.



RECOVERY-M6

- **RECOVERY BED MANUAL**
- 3 Function Bed
- Height/Back & Leg.
- CRCA Perforated Bed Top
- ABS Railings (4Small)
- ABS Head Board
- 125mm Wheel Castors
- Sturdy CRCA Construction.

Options / Add-ons

- **Wheels :** Double Wheel Castors or Central Locking Double Wheel Castors.
- **Bed Top :** 4 Sectional ABS Bed Top & X-Ray Cassette Holder in backrest.
 - **Nurse Control Panel.**
 - **Battery Backup.**
 - **Four Fold Mattress.**



Due to continuous development and innovations, we reserve the rights to change the design and look of the product without prior intimation.